



Clorox Healthcare® Bleach Germicidal Wipes

EPA-registered to kill *C. difficile* spores in 3 minutes
— the fastest contact time available*



Bleach Wipes Most Trusted by U.S. Hospitals

Meet the standards of your healthcare peers and leading institutions:

- 1:10 bleach dilution, 0.55% (5500 ppm) sodium hypochlorite
- Meet CDC and APIC guidelines for surface disinfection of the toughest healthcare pathogens[†]
- Used by 5 of *U.S. News and World Report's* Top 10 Most Influential Hospitals
- Ready to use and premoistened to ensure proper concentration every time
- One-step cleaning and disinfecting; no precleaning required

Fast-acting and effective against 51 microorganisms, including:

- *C. difficile* spores in 3 minutes*
- Representatives of ESKAPE pathogens in 30 seconds, including *Acinetobacter baumannii*, *Enterobacter aerogenes*, VRE and MRSA
- Norovirus in 1 minute
- TB in 3 minutes

* Based on Federal EPA registrations as of July 1, 2012 (3-minute contact time). Use as directed on pre-cleaned, hard, nonporous surfaces.

† Based on Federal master label usage instructions/contact times for *E. coli*, MRSA, *P. aeruginosa*, *S. enterica* and *S. aureus*.

Efficiency and value

Our portfolio of wipes helps your team do the job fast and avoid rework. How?

- **Real-World Contact Times:** Clorox Healthcare® Bleach Germicidal Wipes disinfect quickly and each wipe can keep a 4 sq. ft. area effectively wet for over 2 minutes and 30 seconds:

Clinical Wipe 6" x 5"	Multipurpose Wipe 6.75" x 9"	Terminal Wipe 12" x 12"
2 min 30 sec	3 min 30 sec	5 min

In these cases, no new wipe is required to apply extra wetness.

- **Excellent Surface Coverage:** Clorox Healthcare® Bleach Germicidal Wipes can keep the following surface areas effectively wet for 2 minutes:

Clinical Wipe 6" x 5"	Multipurpose Wipe 6.75" x 9"	Terminal Wipe 12" x 12"
7 sq. ft.	20 sq. ft.	49 sq. ft.

- **Extra-Heavy, Durable Substrate:** Clorox Healthcare® Bleach Germicidal Wipes will not easily tear or fall apart and are designed for better cleaning power. Less tearing means each wipe and each canister lasts longer.



Designed for user and patient comfort

Clorox Healthcare® Bleach Germicidal Wipes are made with a proprietary odor-masking formula for patient and staff comfort. There are several reasons your team can use bleach with confidence.

1. Multiple leading hospital studies noted that when patients and Environmental Services personnel were exposed to bleach, little or no concern about the odor of bleach was expressed!
2. Bleach's odor is caused by the chemical reaction that occurs when bleach begins to break down proteins, such as those that make up infection-causing pathogens.
3. When bleach is used as directed, the chlorine species emitted are lower than recognized levels that cause respiratory irritation and are below the EPA's occupation exposure limits for chlorine.

Broad surface compatibility

To minimize hypochlorite's impact on surfaces, Clorox Healthcare® Bleach Germicidal Wipes are specially formulated with an anti-corrosion agent and are compatible on the following:

- **Surfaces:** stainless steel, plastics, glazed ceramics, glass, porcelain and other materials
- **Areas:** bed rails, tables, equipment surfaces, countertops, floors, toilets, sinks, trash cans, desktop keyboards, telephone receivers, light switches, desks and mobile devices such as IV stands, carts and glucometers

To ensure that the salt left behind by bleach does not build up and cause surface damage over time, wipe down these surfaces after disinfecting with a fresh, clean, damp cloth.

1. Patient and Environmental Service employee satisfaction and tolerance of using germicidal bleach wipes for patient room cleaning. Presentation by Kimberly Aronhalt, BSN, RN, Infection Preventionist, Mayo Clinic, Rochester, Minn., at APIC, New Orleans, July 13, 2010.
 Mayfield JL, Leet T, Miller J, Mundy LM. Environmental Control to Reduce Transmission of *Clostridium difficile*. *Clinical Infectious Diseases*. 2000;31:995-1000.
 Eckstein BC, Adams DA, Eckstein EC et al. Reduction of *Clostridium difficile* and vancomycin-resistant enterococcus contamination of environmental surfaces after an intervention to improve cleaning methods. *BMC Infectious Diseases*. 2007;7:61.

The right-sized wipes for every job

Clorox Healthcare® Bleach Germicidal Wipes come in three different sizes to ensure the right tool for every job. Each size is specifically designed to deliver a perfect clean.

CLINICAL SURFACES	MULTIPURPOSE	PORTABLE	TERMINAL CLEANING
6" x 5"	6.75" x 9"	6.75" x 9"	12" x 12"
			
DESIGNED FOR: Nurses	DESIGNED FOR: Nurses Environmental Services	DESIGNED FOR: Nurses	DESIGNED FOR: Environmental Services
PERFECT FOR: <ul style="list-style-type: none"> Equipment disinfection: Everyday disinfection of clinical surfaces such as glucometers and patient monitors by nurses. 	PERFECT FOR: <ul style="list-style-type: none"> Versatility: Disinfection and cleaning of high-touch surfaces and medical equipment. 	PERFECT FOR: <ul style="list-style-type: none"> Equipment disinfection: Everyday disinfection of clinical surfaces such as glucometers and patient monitors by nurses. 	PERFECT FOR: <ul style="list-style-type: none"> Terminal cleaning
<ul style="list-style-type: none"> High frequency: High consumption and frequent usage. 	<ul style="list-style-type: none"> Spot cleaning and spills 	<ul style="list-style-type: none"> Usage on: Workstations on wheels, transport vehicles, small medical equipment and home health 	<ul style="list-style-type: none"> Reducing labor costs and protecting EVS from spills: Ready-to-use application eliminates formula mixing.
<ul style="list-style-type: none"> Small surfaces: Disinfection and cleaning of high-touch surfaces and medical equipment, such as: Glucometers BP monitors IV pumps Stethoscopes Remote controls Exterior surface of ultrasound transducers Telephone receivers Desktop keyboards Nurse call cords Similar handheld devices 	<ul style="list-style-type: none"> Everyday patient room cleaning: Carts Wheelchairs Countertops Desktops Bed rails Dialysis machines 	<ul style="list-style-type: none"> Soft packs and individually wrapped wipes are ideal for: Infusion pumps Vital signs monitors Portable patient equipment Isolation rooms Workstations on wheels Transport carts 	<ul style="list-style-type: none"> Large surfaces: Cleaning and disinfecting large areas such as: Patient bed area Counters Windowsills Table tops Patient chairs Air vents Toilets Showers Sinks



Use as directed on hard, nonporous surfaces. For any use site that includes electronics, refer to manufacturers' instructions before use.

Clorox Healthcare® Bleach Germicidal Wipes

Effective against: Kill Time

	Effective against:	Kill Time
Bacteria	<i>Acinetobacter baumannii</i>	30 sec
	<i>Bordetella pertussis</i>	30 sec
	<i>Campylobacter jejuni</i>	30 sec
	<i>Clostridium difficile</i> spores	3 min
	<i>Enterobacter aerogenes</i>	30 sec
	<i>Enterobacter cloacae</i>	30 sec
	Vancomycin-resistant <i>Enterococcus faecalis</i> (VRE)	30 sec
	Multidrug-resistant <i>Enterococcus faecium</i>	30 sec
	<i>Escherichia coli</i>	30 sec
	<i>Escherichia coli</i> O157:H7	30 sec
	ESBL-producing <i>Escherichia coli</i>	30 sec
	<i>Escherichia coli</i> NDM-1	30 sec
	<i>Klebsiella pneumoniae</i>	30 sec
	ESBL-producing <i>Klebsiella pneumoniae</i>	30 sec
	Carbapenem-resistant <i>Klebsiella pneumoniae</i>	30 sec
	<i>Klebsiella pneumoniae</i> NDM-1	30 sec
	<i>Legionella pneumophila</i>	30 sec
	<i>Listeria monocytogenes</i>	30 sec
	<i>Proteus mirabilis</i>	30 sec
	<i>Pseudomonas aeruginosa</i>	30 sec
	<i>Salmonella enterica</i>	30 sec
	<i>Serratia marcescens</i>	30 sec
	<i>Shigella dysenteriae</i>	30 sec
	<i>Staphylococcus aureus</i>	30 sec
Linezolid-resistant <i>Staphylococcus aureus</i> (LRSA)	30 sec	
Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	30 sec	
Methicillin-resistant <i>Staphylococcus aureus</i> (CA-MRSA) (USA 300)	30 sec	
Methicillin-resistant <i>Staphylococcus aureus</i> (CA-MRSA) (USA 400)	30 sec	
Vancomycin-resistant <i>Staphylococcus aureus</i> (VRSA)	30 sec	
Vancomycin intermediate-resistant <i>Staphylococcus aureus</i> (VISA)	30 sec	
<i>Streptococcus pneumoniae</i>	30 sec	
<i>Streptococcus pyogenes</i>	30 sec	
Viruses	Adenovirus Type 2	1 min
	Avian Influenza A Virus	1 min
	Hepatitis A Virus (HAV)	1 min
	Hepatitis B Virus (HBV)	1 min
	Hepatitis C Virus (HCV)	1 min
	Herpes Simplex Virus (HSV-2)	1 min
	Human Immunodeficiency Virus Type 1 (HIV-1)	30 sec
	Human Coronavirus	1 min
	Influenza A Virus	1 min
	Norovirus*	1 min
	Poliovirus	1 min
	Respiratory Syncytial Virus (RSV)	1 min
	Rhinovirus	1 min
Rotavirus	1 min	
Fungi, Vet	<i>Candida albicans</i>	3 min
	Canine Parvovirus	3 min
	Feline Panleukopenia Virus	3 min
TB	<i>Trichophyton mentagrophytes</i> (Athlete's Foot Fungus)	3 min
	<i>Mycobacterium bovis</i>	3 min

TOTAL: 51 Claims

Use as directed on hard, nonporous surfaces.
 * Obtained via surrogate using Feline calicivirus.

For more information, contact your Clorox sales representative or call **800-234-7700**.
 email: healthcare@clorox.com
 visit us: www.cloroxhealthcare.com

© 2015 Clorox Professional Products Company,
 1221 Broadway, Oakland, CA 94612.

Clorox Healthcare® Bleach Germicidal Wipes EPA Reg. No. 67619-12

Clinical Wipes

6" x 5" Wipes
 6/150 ct.
 Item No. 30577



Multipurpose Wipes

6.75" x 9" Wipes
 6/70 ct.
 Item No. 35309



Terminal Wipes

12" x 12" Wipes
 2/110 ct. Bucket
 Item No. 30358



Terminal Wipes Refill

12" x 12" Wipes
 2/110 ct. Refill
 Item No. 30359



Germicidal Wipes

6.75" x 9" Wipes
 Transport Pack 24/20 ct.
 Item No. 31469
 Individual 6/50 ct.
 Item No. 31424

